

## Manufacturers Instrument Compatibility

CBI 1501 - Creative Biomedics – 1702 WinDX  
Alliance – InCheck # 3109750  
Alliance - InCheck Dial - # 3109912  
Micro-Direct – MicroLab # ML3500S  
Micro-Direct – Spiro - # ML2525  
Micro-Direct – Micro # MS01  
Micro-Direct - Micro Plus # MS03  
Micro-Direct – Micro Plus # MS04  
Micro-Direct – MicroGP # MS07  
Micro-Direct – Micro Diary # MD02  
Micro-Direct – Micro Loop # ML3535  
Vacu Med – Discovery 2 - # 16503  
Vacu Med – Discovery 2 - # 16504  
Vacu Med – Spirolab - #16200  
Vacu Med - Spirolab 2 - #16210  
Vacu Med – Wizard – #S1000  
Vacu Med – Micro Loop - #16700  
Vacu Med – Micro GP - #S-1500  
Vacu Med – Micro DL - # S-1575  
Vacu Med – Micro – # S-1200  
Vacu Med – Micro Plus - # S-1275  
CDX Corp – SpiroLab 2 - #  
Future Med – SpiroPerfect- # WASPIROS  
Welch Allyn – SpiroPerfect - # SPIROS  
Vitalograph – Vitalograph 2120  
Vitalograph – Vitalograph – Vitalograph Alpha  
Vitalograph – Pneumotrac  
SDI Diagnostics - Spirolab 2 - # **29-1015**  
**SDI Diagnostics - Spirolab - # 29-1010**  
**SDI Diagnostics – SPG Muti Functions - 29-1003**  
**Shiller – SP 1 - #954000**  
**MIR – Spirobank #**

CBI – 1502 – Alliance - All Flow #550100  
Alliance – One Flow #3141001  
Viasys – VMAX Encore  
Viasys – SensorMedics #2200  
Viasys – SensorMedics - #2400  
Viasys – Sensor Medics - # 2450  
Viasys – Sensor Medics - #2800

CBI-1503 nSpire Health – KoKo #313100  
KoKo #313105  
KoKo Legend - #314000  
KoKo Mate #355100  
KoKo Trek #333300  
KoKo Trek #333305  
Eaglet - #K003510  
**Assess** Peak Flow Meter

**Burdick** – Spirosense

**Circadian** – All Models (if using black Tamarac pneumotach, need 7975)

**Cybermedics** – Spinnaker TL, SL, CM3, CM310 (if white or black cone use 7980)

**Cybermedics** – Moose

**Godart** – All Models

**Gould** – 5004

**Healthscan** – Assess Peak Flow Meter

**Jaeger** – All Models except Spiro Pro

**Jones** – All Models except Satellite (small end of filter)

**Koko** – All Models

**Med Graphics** – 1070, 1085 (w/Prevent, use PG II; if no Prevent, 7970-spirometry, 7961-vol/diff)

**Med Science** – 570, 6257

**Med Science** – 3000, 4000 (spirometry-no adapt./vol-7961/diff-7960B)

**Med Systems** – Spirolite 101 & 201 (small end of filter), 343

**New Med** – 7000 I (Spirobyte)

**PK Morgan** – Benchmark (spirometry-7978, vol&diff-7960B), Flexi Flow (small end of filter), Trans Flow

**Pneumedics** – All Models (small end of filter)

**Riko** – AS300

**S&M** – CF-150B, GS520, Keystone III, Portascreen, Pulmoscreen

**ViaSys/Sensormedics** – Horizon, Ohio, 922, 2130, 2200, 2400, 2450, 2800, 6200, 2100 (spirometry-7961, vol-7960B, diff-7960B), Vmax (can use PGII w/7976 or 7977)

**Spirolink** – All Models

**Spirometrics** – LTE, 2500 LTE

**Spiroscreen** – All Models (small end of filter)

**Spirotech** – All Models

**Tamarac** – All